

Glucagon ELISA Kit For Mouse, Rat & Human Samples

**10 μ L
Sample**

■ **HIGHLY
SPECIFIC**

■ **USES
10 μ L
SAMPLE**

■ **GREAT
PRICE!**

Catalog Number

Mouse Glucagon ELISA Kit81518
Rat Glucagon ELISA Kit81519
Human Glucagon ELISA Kit81520

Sample Size

10 μ L

Sample Type

Plasma, Serum, Cell Culture

Sensitivity

< 1.1 pg/mL
(0.3 pmol/L)

Range

7.8 - 500 pg/mL
(2.2 - 144 pmol/L)

No. Tests

96 Tests

Highlights

- Small sample volume (10 μ L)
- High sensitivity (< 1.1 pg/mL or 0.3 pmol/L)
- 1-step ELISA
- Short, overnight procedure
- Under 1% cross-reactivity with glicentin, oxyntomodulin, mini-glucagon, GIP, GLP-1, or GLP-2 (7-36 or 9-36)NH₂

Order online at www.crystalchem.com