

Mouse and Rat Creatinine Kits

Lower Interference vs Jaffe Method

Catalog Numbers

Mouse Creatinine Kit..... 80350

Rat Creatinine Kit..... 80340

Sample Size	Assay time	Precision
8 μ L	< 15 mins	CV < 10%

Samples	Range
Serum, Urine	0.15 - 13.5 mg/dL (Serum)
	0.15 - 141 mg/dL (Urine)

Highlights

- Enzymatic method results in lower interference vs Jaffe method
- Uses only 8 μ L sample volume
- Fast (< 15 minutes)

Order Online
www.crystalchem.com