

# Mouse and Rat Glucose Assay Kits

**Fast  
Assay!**

■ **SMALL  
SAMPLE  
VOLUME**

■ **FAST**

■ **HIGHLY  
ACCURATE**

## Catalog Numbers

Mouse Glucose Assay Kit.....	81692
Rat Glucose Assay Kit.....	81693

Sample Size	Range	Precision
2 $\mu$ L	0 - 675 mg/dL	CV < 10%

Sample Type	No. Tests	Assay Time
Plasma, Serum,	96 Tests	5 mins
Urine		

## Highlights

- Very fast (5 min)
- Uses only 2  $\mu$ L sample volume
- Can be used with serum, plasma, and urine
- More accurate results than glucometer for mouse and rat serum/plasma
- Human kit (#81696) and 10-packs are also available

Order online at [www.crystalchem.com](http://www.crystalchem.com)