

Ultra Sensitive Mouse Glucagon ELISA Kit

**10 μ L
Sample**

- **10 μ L Sample**
- **Easy-to-use**
- **High Sensitivity**
- **1-Step ELISA**
- **Great price!**

Catalog Number

NEW Mouse Glucagon ELISA Kit81518

Sample Size

10 μ L

Sample Type

Plasma, Serum, Cell Culture

Sensitivity

< 1.1 pg/mL
(0.3 pmol/L)

Range

7.8 - 500 pg/mL
(2.2 - 144 pmol/L)

No. Tests

96 Tests

Highlights

- Small sample volume (10 μ L)
- High sensitivity (< 1.1 pg/mL or 0.3 pmol/L)
- 1-step ELISA
- Short, overnight procedure
- Under 1% cross-reactivity with glicentin, oxyntomodulin, mini-glucagon, GIP, GLP-1, or GLP-2 (7-36 or 9-36)NH₂
- Also available for Rat and Human samples

Order online at www.crystalchem.com