

Ultra Sensitive Mouse Insulin ELISA Kit, 2nd Gen

**NEW
KIT!**

Liquid standards

5 μ L sample

One wash step

No cross reactivity to:

-> Mouse proinsulin

-> Human insulin

Catalog Number

Ultra Sensitive Mouse Insulin
ELISA Kit, 2nd Gen.....Cat# 62100

Sample Size	Range	Precision
5 μ L	0.09 - 11 ng/mL	CV \leq 10%

Samples	Sensitivity
Serum, plasma, cell culture	\leq 0.05 ng/mL

Highlights

- Easy to use with all liquid standards
- Highly sensitive with only 5 μ L of sample
- No cross reactivity to mouse proinsulin or human insulin
- One wash step
- Wide dynamic range of 0.09 - 11 ng/mL
- 10-kit bundle available (Cat# 62105)

Order online at www.crystalchem.com