

Ultra Sensitive Rat Insulin ELISA Kit, 2nd Gen

**NEW
KIT!**

Liquid standards

5 μ L sample

One wash step

No cross reactivity to:

-> Rat proinsulin

-> Human insulin

Catalog Number

Ultra Sensitive Rat Insulin

ELISA Kit, 2nd Gen.....Cat# 62300

Sample Size

5 μ L

Range

0.1 - 13 ng/mL

Precision

CV \leq 10%

Samples

Serum, plasma, cell culture

Sensitivity

\leq 0.05 ng/mL

Highlights

- Easy to use with all liquid standards
- Highly sensitive with only 5 μ L of sample
- No cross reactivity to rat proinsulin or human insulin
- One wash step
- Wide dynamic range of 0.1 - 13 ng/mL
- 10-kit bundle available (Cat# 62305)

Order online at www.crystalchem.com